

4-Formylbenzoic acid


Catalog No: #S1691

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Formylbenzoic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	150.13
Formulation	C ₈ H ₆ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)C1=CC=C(C=O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Xavier TS, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2013 Oct;114:502-8.

Note: This product is for in vitro research use only