

## 7-Hydroxycoumarin


Catalog No: #S1693

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	7-Hydroxycoumarin
Brief Description	Inhibitors
Purification	98%+
Target Name	Antioxidant
Calculated MW	162.14
Formulation	C <sub>9</sub> H <sub>6</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: OC1=CC2=C(C=C1)C=CC(=O)O2  
Pathways: oxidation-reduction  
Receptor: Antioxidant  
Boiling pt: Melting  
pt: Solubility: Soluble in DMSO  
Appearance: Gray Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Kabeya LM, et al. Chem Biol Interact. 2013 Oct 25;206(1):63-75.

Note: This product is for in vitro research use only