


Description

Product Name	DL-O-Tyrosine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives Agonist
Calculated MW	181.19
Formulation	C9H11NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:NC(CC1=CC=CC=C1O)C(O)=OPathways:OthersReceptor:Amino Acids and DerivativesBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Koyama H. J Biochem. 1984 Aug;96(2):421-7.

Note: This product is for in vitro research use only