

debutylbupivacaine

Catalog No: #S1697


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	debutylbupivacaine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	232.32
Formulation	C14H20N2O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC1=CC=CC(C)=C1NC(=O)C1CCCN1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pihlajamäki K, et al. Br J Anaesth. 1990 May;64(5):556-62.

Note: This product is for in vitro research use only