

Carvedilol EP Impurity E

Catalog No: #S1698


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carvedilol EP Impurity E
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	203.67
Formulation	C ₉ H ₁₃ NO ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: Cl.COC1=CC=CC=C1OCCN
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only