

## Monoethylglycinexylidide

Catalog No: #S1700


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Monoethylglycinexylidide
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	206.28
Formulation	C <sub>12</sub> H <sub>18</sub> N <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Other SMILES:CCNCC(=O)NC1=C(C)C=CC=C1C Pathways:Others Receptor:Others Boiling pt:Melting pt:Solubility:Limited solubility Appearance:White Powder Remark:For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Gerlach JC, et al. Artif Organs. 2010 Jun;34(6):462-72.

Note: This product is for in vitro research use only