


Description

Product Name	H-ASP-PHE-OH
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	280.28
Formulation	C13H16N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1=CC=C(C=C1)C[C@@H](C(=O)O)NC(=O)[C@H](CC(=O)O)N
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Limited solubility
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Murata T, et al. J Biotechnol. 1993 Apr;28(2-3):301-12.

Note: This product is for in vitro research use only