

Nortropinone hydrochloride

Catalog No: #S1702


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Nortropinone hydrochloride |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 161.63 |
| Formulation | C7H12ClNO |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
 SMILES: Cl.O=C1CC2CCC(C1)N2
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Limited
 solubility
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Tea I, et al. Anal Bioanal Chem. 2011 Sep; 401(4):1263-71.

Note: This product is for in vitro research use only