

4-Amino-5-imidazolecarboxamide


Catalog No: #S1705

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Amino-5-imidazolecarboxamide
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	126.12
Formulation	C ₄ H ₆ N ₄ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: NC(=O)C1=C(N)N=CN1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited
solubility Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pastore S, et al. Rheumatol Int. 2015 Apr;35(4):619-27.

Note: This product is for in vitro research use only