

## 4-Hydroxy-3-nitrophenylacetic acid


Catalog No: #S1706

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	4-Hydroxy-3-nitrophenylacetic acid
Brief Description	Inhibitors
Purification	96.73%
Target Name	Others
Calculated MW	197.14
Formulation	C <sub>8</sub> H <sub>7</sub> NO <sub>5</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: OC(=O)CC1=CC(=C(O)C=C1)[N+](=O)[O-]  
Pathways: Others  
Receptor: Ig gamma-2 chain C region; Ig heavy chain V-I region  
Boiling pt: Melting pt: Solubility: Limited solubility  
Appearance: Yellow Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Berman HM, et al. Nucleic Acids Res. 2000 Jan 1;28(1):235-42.

Note: This product is for in vitro research use only