

2'-Deoxyinosine


Catalog No: #S1709

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2'-Deoxyinosine
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	252.23
Formulation	C10H12N4O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system SMILES: C1[C@@H]([C@H](O[C@H]1N2C=NC3=C2NC(=O)C(=O)N3)CO)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Baumann T, et al. BMC Biotechnol. 2013 Oct 4;13:81.

Note: This product is for in vitro research use only