


Description

Product Name	Echinacoside
Brief Description	Inhibitors
Purification	100.00%
Target Name	Others
Calculated MW	786.73
Formulation	C35H46O20
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@H]2[C@H]([C@@H](O[C@H]([C@H]2OC(=O)/C=C/C3=CC(=C(C=C3)O)O)CO[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)OC(=O)C5=CC(=C(C=C5)O)O)O)OP

Pathways:OthersReceptor:OthersBoiling
pt:Melting pt:Solubility:DMSO: 10 mMAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhang D, et al. Int J Biol Macromol. 2014 Sep 2;72C:243-253.

Note: This product is for in vitro research use only