

SCR7 pyrazine

Catalog No: #S1724

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SCR7 pyrazine
Brief Description	Inhibitors
Purification	99.85%
Target Name	DNA/RNA Synthesis inhibitor
Calculated MW	332.38
Formulation	C18H12N4OS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: O=C1NC(=S)NC2=C1N=C(C1=CC=CC=C1)C(=N2)C1=CC=CC=C1 Pathways: DNA Damage/DNA

Repair Receptor: DNALigase IV Boling pt: Melting pt: Solubility: DMSO: 66 mg/mL (198.56 mM) Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Srivastava M, et al. Cell. 2012, 151(7), 1474-1487.

Note: This product is for in vitro research use only