

TOS-ARG-OME hydrochloride

Catalog No: #S1731


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TOS-ARG-OME hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives inhibitor
Calculated MW	378.87
Formulation	C ₁₄ H ₂₃ CIN ₄ O ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: CC1=CC=C(C=C1)S(=O)(=O)N[C@@H](CCCN=C(N)N)C(=O)OC
 Pathways: Others
 Receptor: Amino Acids and Derivatives
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only