

Rabeprazole Related Compound E

Catalog No: #S1734


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabeprazole Related Compound E
Brief Description	Inhibitors
Purification	99.41%
Target Name	Proton pump inhibitor
Calculated MW	343.44
Formulation	C18H21N3O2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:COCCOC1=C(C)C(CSC2=NC3=CC=CC=C3N2)=NC=C1Pathways:Membrane transporter/Ion channelReceptor:Proton pumpBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ohara T, et al. Helicobacter. 2001 Jun;6(2):125-9.

Note: This product is for in vitro research use only