

## Apixaban

Catalog No: #S1736


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Apixaban
Brief Description	Inhibitors
Purification	99.54%
Target Name	Factor Xa inhibitor
Calculated MW	459.5
Formulation	C25H25N5O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism

systemSMILES:COC1=CC=C(C=C1)N1N=C(C(N)=O)C2=C1C(=O)N(CC2)C1=CC=C(C=C1)N1CCCCC1=OPathways:MetabolismReceptor:Factor Xa (human); Factor Xa (rabbit)Boiling pt:Melting pt:Solubility:DMSO: 18 mg/mL (39.17 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Pinto DJ, et al. J Med Chem. 2007, 50(22), 5339-5356.

Note: This product is for in vitro research use only