


Description

Product Name	WZ811
Brief Description	Inhibitors
Purification	99.93%
Target Name	CXCR antagonist
Calculated MW	290.36
Formulation	C18H18N4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=CC=C(NC2=CC=CC=N2)C(=C(C=C1)C3=CC=CC=N3)
Pathways: Autophagy; GPCR/G Protein;
Immunology/Inflammation Receptor: CXCR4
Boiling pt: Melting pt: Solubility: DMSO: 30 mg/mL (103.32 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhan W, et al. J Med Chem, 2007, 50(23), 5655-5664.

Note: This product is for in vitro research use only