

TAK-063


Catalog No: #S1745

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TAK-063
Brief Description	Inhibitors
Purification	99.00%
Target Name	PDE inhibitor
Calculated MW	428.42
Formulation	C23H17FN6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
SMILES: COC1=CN(N=C(C2=CC=NN2C2=CC=CC=C2)C1=O)C1=C(F)C=C(C=C1)N1C=CC=N1
Pathways: Metabolism
Receptor: PDE10A
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kunitomo J, et al. J Med Chem. 2014 Nov 26;57(22):9627-43.

Note: This product is for in vitro research use only