

Torezolid


Catalog No: #S1747

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Torezolid
Brief Description	Inhibitors
Purification	96.97%
Target Name	MAO inhibitor
Calculated MW	370.34
Formulation	C17H15FN6O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CN1N=C(N=N1)C2=NC=C(C=C2)C3=C(C=C(C=C3)N4C[C@@H](OC4=O)CO)F Pathways: Metabolism Receptor: MAO-A; MAO-B
Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Solid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Choi S, et al. Antimicrob Agents Chem Others. 2012 Jun 19.

Note: This product is for in vitro research use only