


Description

Product Name	AL082D06
Brief Description	Inhibitors
Purification	98.00%
Target Name	GR antagonist
Calculated MW	409.91
Formulation	C23H24ClN3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CN(C)C1=CC=C(C=C1)C(C1=CC=C(C=C1)N(C)C)C1=C(Cl)C=CC(=C1)[N+](O-)=O

Pathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Avenant C, et al. Mol Cell Endocrinol. 2010 Oct 7;327(1-2):72-88.

Note: This product is for in vitro research use only