

Eliprodil

Catalog No: #S1751


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eliprodil
Brief Description	Inhibitors
Purification	99.14%
Target Name	NMDAR antagonist
Calculated MW	347.85
Formulation	C ₂₀ H ₂₃ ClFNO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:OC(CN1CCC(CC2=CC=C(F)C=C2)CC1)C1=CC=C(Cl)C=C1Pathways:NeuroscienceReceptor:NR2B-NMDABoling pt:Melting pt:Solubility:DMSO: 25 mMAppearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Bath CP, et al. Eur J Pharmacol. 1996 Mar 28;299(1-3):103-12.

Note: This product is for in vitro research use only