

## PD153035 hydrochloride

Catalog No: #S1761


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PD153035 hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	EGFR inhibitor
Calculated MW	396.67
Formulation	C16H15BrClN3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: Cl.COC1=C(OC)C=C2C(NC3=CC(Br)=CC=C3)=NC=NC2=C1  
Pathways: Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/Adaptors Receptor: EGFR; EGFR  
Boiling pt: Melting pt: Solubility: DMSO: 0.5 mg/mL (1.26 mM)  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Fry DW, et al. Science. 1994, 265(5175), 1093-1095.

Note: This product is for in vitro research use only