

Ro 48-8071 fumarate

Catalog No: #S1771


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ro 48-8071 fumarate
Brief Description	Inhibitors
Purification	98.72%
Target Name	Others inhibitor
Calculated MW	564.44
Formulation	C27H31BrFNO6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: CN(CCCCCOC1=CC(=C(C=C1)C(=O)C2=CC=C(C=C2)Br)F)CC=C.C(=C/C(=O)O)\C(=O)O
 Pathways: Others
 Receptor: Oxidosqu
 alene cyclase
 Boling pt:
 Melting pt:
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Morand OH, et al. J Lipid Res. 1997 Feb;38(2):373-90.

Note: This product is for in vitro research use only