

Apoptosis Activator 2

Catalog No: #S1772


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apoptosis Activator 2
Brief Description	Inhibitors
Purification	98.00%
Target Name	Caspase activator
Calculated MW	306.14
Formulation	C15H9Cl2NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C1C=C(Cl)C=C(CN2C(=O)C(=O)C3=CC=CC=C23)C=C1Pathways:Apoptosis;

Proteases/ProteasomeReceptor:Caspase-3Boiling pt:Melting pt:Solubility:DMSO: 61 mg/mL (199.25 mM)Appearance:Orange PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nguyen JT, et al. Proc Natl Acad Sci U S A, 2003, 100(13), 7533-7538.

Note: This product is for in vitro research use only