


Description

Product Name	LXR623
Brief Description	Inhibitors
Purification	99.00%
Target Name	Liver X Receptor agonist
Calculated MW	422.78
Formulation	C21H12ClF5N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: FC1=CC=C(C=C1)C1=C2C=CC=C(C2=NN1CC1=C(Cl)C=C(F)C=C1)C(F)(F)F
 Pathways: Others
 Receptor: LXR α ; LXR β
 Boiling pt: Melting pt: Solubility: DMSO: > 10mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Quinet EM, et al. J Lipid Res. 2009 Dec;50(12):2358-70.

Note: This product is for in vitro research use only