

Levobupivacaine

Catalog No: #S1787


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Levobupivacaine
Brief Description	Inhibitors
Purification	99.70%
Target Name	Sodium Channel inhibitor
Calculated MW	288.43
Formulation	C18H28N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCCCN1CCCC[C@H]1C(=O)NC2=C(C=CC=C2C)C
 Pathways: Membrane transporter/Ion channel
 Receptor: Sodium Channel
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ueta K, et al. Reg Anesth Pain Med. 2006 Jan-Feb;31(1):19-25.

Note: This product is for in vitro research use only