

Factor XIII B Polyclonal Antibody

Catalog No: #40911

Package Size: #40911-1 50ul #40911-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Factor XIII B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Factor XIII B Polyclonal Antibody detects endogenous levels of Factor XIII B protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Factor XIII B.
Target Name	Factor XIII B
Other Names	F13B; Coagulation factor XIII B chain; Fibrin-stabilizing factor B subunit; Protein-glutamine gamma-glutamyltransferase B chain; Transglutaminase B chain
Accession No.	Swiss-Prot: P05160NCBI Gene ID: 2165
Uniprot	P05160
GeneID	2165;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/20000

Images



Western Blot analysis of COLO cells using Factor XIII B Polyclonal Antibody

Note: This product is for in vitro research use only