

## HPGDS-inhibitor-1

Catalog No: #S1804


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	HPGDS-inhibitor-1
Brief Description	Inhibitors
Purification	98.00%
Target Name	PGE Synthase inhibitor
Calculated MW	381.37
Formulation	C19H19F4N3O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: FC1=CC=CC(=C1)C1=NC=C(C=C1)C(=O)NC1CCN(CC(F)(F)F)CC1 Pathways: Immunology/Inflammation Receptor: hPGDS

Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Chris P. Carron, et al. ACS Med. Chem. Lett., 2010, 1 (2), pp 59-63

Note: This product is for in vitro research use only