

## FADD Polyclonal Antibody

Catalog No: #40912

Package Size: #40912-1 50ul #40912-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FADD Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	FADD Polyclonal Antibody detects endogenous levels of FADD protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human FADD around the non-phosphorylation site of S194.
Target Name	FADD
Other Names	FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein; FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity
Accession No.	Swiss-Prot: Q13158NCBI Gene ID: 8772
Uniprot	Q13158
GeneID	8772;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

## Images

Western Blot analysis of HeLa cells using FADD Polyclonal Antibody



Note: This product is for in vitro research use only