

Purmorphamine

Catalog No: #S1810


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purmorphamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Hedgehog/Smoothened antagonist
Calculated MW	520.62
Formulation	C31H32N6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C1CCC(CC1)N1C=NC2=C1N=C(OC1=CC=CC3=CC=CC=C13)N=C2NC1=CC=C(C=C1)N1CCOCC1 Pathways: GPCR/G Protein; Stem Cells Receptor: Smoothened Bolding pt: Melting pt: Solubility: DMSO: 4 mg/mL warmed (7.68 mM) Appearance: Faint Green Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Faghihi F, Biomed Pharmacol, 2012, 3322(12).

Note: This product is for in vitro research use only