


## Description

Product Name	KX1-004
Brief Description	Inhibitors
Purification	98.00%
Target Name	Src inhibitor
Calculated MW	284.29
Formulation	C16H13FN2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:OC1=CC=CC(CNC(=O)C2=CC3=C(N2)C=CC(F)=C3)=C1Pathways:AngiogenesisReceptor:Src-PTKBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Harris KC, et al. Hear Res. 2005 Oct;208(1-2):14-25.

Note: This product is for in vitro research use only