

ICILIN

Catalog No: #S1814


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ICILIN
Brief Description	Inhibitors
Purification	98.00%
Target Name	TRP/TRPV Channel agonist
Calculated MW	311.29
Formulation	C16H13N3O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: OC1=CC=CC=C1N1CC=C(NC1=O)C1=CC(=CC=C1)[N+](=[O-])=O
 Pathways: Membrane transporter/Ion channel
 Receptor: TRPM8
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lee S, et al. Biochem Biophys Res Commun. 2013 Nov 29;441(4):1005-10.

Note: This product is for in vitro research use only