

Clemizole

Catalog No: #S1822


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Clemizole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histamine Receptor inhibitor
Calculated MW	325.84
Formulation	C19H20ClN3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: C1C=CC(CN2C(CN3CCCC3)=NC3=CC=CC=C23)C=C1
 Pathways: GPCR/G Protein;
 Immunology/Inflammation; Neuroscience
 Receptor: NS4BH1 histamine receptor
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Einav S, et al. J Infect Dis. 2010 Jul 1;202(1):65-74.

Note: This product is for in vitro research use only