

StemRegenin 1 (SR1)

Catalog No: #S1831


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	StemRegenin 1 (SR1)
Brief Description	Inhibitors
Purification	98.00%
Target Name	AhR antagonist
Calculated MW	429.54
Formulation	C24H23N5OS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC(C)N1C=NC2=C1N=C(N=C2NCCC1=CC=C(O)C=C1)C1=CSC2=CC=CC=C12 Pathways: Metabolism Receptor: AhR Boling
 pt: Melting pt: Solubility: DMSO: 86 mg/mL (200.21 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Boitano AE, et al. Science, 2010, 329(5997), 1345-1348.

Note: This product is for in vitro research use only