

## Walrycin B


Catalog No: #S1832

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Walrycin B
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibacterial inhibitor
Calculated MW	337.26
Formulation	C <sub>14</sub> H <sub>10</sub> F <sub>3</sub> N <sub>5</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CN1N=C(N=C2C(=O)N(C)C(=O)N=C12)C1=CC=C(C=C1)C(F)(F)F Pathways: Microbiology & virology Receptor: bacterial WalR response regulator  
 Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Gotoh Y, et al. J Antibiot (Tokyo). 2010 Mar;63(3):127-34.

Note: This product is for in vitro research use only