


Description

Product Name	UNC0379
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histone Methyltransferase inhibitor
Calculated MW	413.56
Formulation	C23H35N5O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:COC1=C(OC)C=C2C(NCCCCCN3CCCC3)=NC(=NC2=C1)N1CCCC1Pathways:Chromatin/EpigeneticReceptor:SETD8Boling

pt:Melting pt:Solubility:DMSO: > 10mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ma A, et al. J Med Chem. 2014, 57(15), 6822-6833.

Note: This product is for in vitro research use only