

B-Raf IN 1

Catalog No: #S1845


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	B-Raf IN 1
Brief Description	Inhibitors
Purification	98.00%
Target Name	Raf inhibitor
Calculated MW	515.53
Formulation	C29H24F3N5O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CN(C)CC1=CC=C(C=C1)C1=C2N=CC=C(N2N=C1)C1=CC(NC(=O)C2=CC(=CC=C2)C(F)(F)F)=CC=C1 Pathways: MAPK Receptor
 :B-Raf; C-Raf Boiling pt: Melting pt: Solubility: DMSO: 40 mg/mL Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Berger DM, et al. Bioorg Med Chem Lett. 2009 Dec 1;19(23):6519-23.

Note: This product is for in vitro research use only