

Necrostatin-1


Catalog No: #S1847

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Necrostatin-1
Brief Description	Inhibitors
Purification	99.25%
Target Name	TNF inhibitor
Calculated MW	259.33
Formulation	C13H13N3OS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CN1C(=S)NC(CC2=CNC3=CC=CC=C23)C1=O
Pathways: Apoptosis Receptor: RIP1
Boiling pt: Melting pt: Solubility: DMSO: 51 mg/mL warmed (196.66 mM)
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Han W, et al. Apoptosis. 2009 May;14(5):674-86.

Note: This product is for in vitro research use only