


Description

Product Name	TPT-260
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	333.32
Formulation	C ₈ H ₁₄ Cl ₂ N ₄ S ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Other SMILES:Cl.Cl.NC(=N)SCC1=CC=C(CSC(N)=N)S1 Pathways:Others Receptor:Others Boiling pt:Melting pt:Solubility:DMSO: 10 mM Appearance:White Powder Remark:For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only