


Description

Product Name	CID 797718
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	217.22
Formulation	C12H11NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC1=CC2=C(OC(=O)C3=C2CCCN3)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only