


## Description

Product Name	AT-56
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others inhibitor
Calculated MW	397.52
Formulation	C25H27N5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:C(CCC1=NNN=N1)CN1CCC(CC1)=C1C2=CC=CC=C2C=CC2=CC=CC=C12Pathways:OthersReceptor:L-PGDSBoiling pt:Melting pt:Solubility:DMSO: > 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Urbanet R, et al. Hypertension. 2015 Jul;66(1):149-57.

Note: This product is for in vitro research use only