

Fenretinide

Catalog No: #S1872


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fenretinide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Retinoid Receptor agonist
Calculated MW	391.55
Formulation	C26H33NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CC1=C(C(CCC1)(C)C)/C=C/C(=C/C=C/C(=C/C(=O)NC2=CC=C(C=C2)O)/C)/CPathways:MetabolismReceptor:RAR/RAXBoling
pt:Melting pt:Solubility:DMSO: 100mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Manzano VM, et al. J Pharmacol Exp Ther, 1999, 289(1), 123-132.

Note: This product is for in vitro research use only