

FAM80B Polyclonal Antibody

Catalog No: #40919

Package Size: #40919-1 50ul #40919-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FAM80B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt Mk
Specificity	FAM80B Polyclonal Antibody detects endogenous levels of FAM80B protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human FAM80B.
Target Name	FAM80B
Other Names	RIMKLB; FAM80B; KIAA1238; Beta-citryl-glutamate synthase B; N-acetyl-aspartyl-glutamate synthetase B; NAAG synthetase B; NAAGS; Ribosomal protein S6 modification-like protein B
Accession No.	Swiss-Prot: Q9ULI2NCBI Gene ID: 57494
Uniprot	Q9ULI2
GeneID	57494;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HepG2 cells using FAM80B Polyclonal Antibody

Note: This product is for in vitro research use only