

Atglistatin


Catalog No: #S1875

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atglistatin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Lipase inhibitor
Calculated MW	283.37
Formulation	C17H21N3O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: CN(C)C(=O)NC1=CC=CC(=C1)C1=CC=C(C=C1)N(C)C
Pathways: Metabolism
Receptor: ATGL
Boiling pt: Melting pt: Solubility: Ethanol: 4 mg/mL (14.11 mM); DMSO: 57 mg/mL (201.15 mM)
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Mayer N, et al. Nat Chem Biol. 2013, 9(12), 785-787.

Note: This product is for in vitro research use only