

Methylisoidigotin

Catalog No: #S1882


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methylisoidigotin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Apoptosis inducer inhibitor
Calculated MW	276.29
Formulation	C17H12N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CN1C2=CC=CC=C2/C(=C\C3/C4=CC=CC=C4NC3=O)/C1=OPathways:ApoptosisReceptor:Apoptosis inducerBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Dark Purple PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lee CC, et al. Leuk Lymphoma. 2010 May;51(5):897-905.

Note: This product is for in vitro research use only