


Description

Product Name	Inauhzin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sirtuin inhibitor
Calculated MW	469.58
Formulation	C25H19N5OS2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CCC(SC1=NC2=C(N=N1)C1=CC=CC=C1N2)C(=O)N1C2=CC=CC=C2SC2=CC=CC=C12Pathways:Chromatin/Epigenetic; DNA Damage/DNA RepairReceptor:SIRT1Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhang Q, et al. EMBO Mol Med. 2012 Apr; 4(4): 298-31

Note: This product is for in vitro research use only