


## Description

Product Name	CCT007093
Brief Description	Inhibitors
Purification	98.00%
Target Name	phosphatase inhibitor
Calculated MW	272.39
Formulation	C15H12OS2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:C1/C(=C/C2=CC=CS2)/C(=O)/C(=C/C3=CC=CS3)/C1Pathways:OthersReceptor:PPM1DBoling pt:Melting pt:Solubility:DMSO: 1 mg/mL warmed (3.67 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Lee JS, et al. J Dermatol Sci. 2014, 73(2), 125-134.

Note: This product is for in vitro research use only