

Filgotinib (GLPG0634)

Catalog No: #S1929


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Filgotinib (GLPG0634)
Brief Description	Inhibitors
Purification	99.00%
Target Name	JAK inhibitor; Tyrosine Kinases inhibitor
Calculated MW	425.5
Formulation	C21H23N5O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: O=C(NC1=NN2C(C=CC=C2C2=CC=C(CN3CCS(=O)(=O)CC3)C=C2)=N1)C1CC1
 Pathways: Angiogenesis; Proteases/Proteasome; Chromatin/Epigenetic; JAK/STAT signaling; Stem Cells
 Receptor: JAK1; JAK2; JAK3; Tyk2
 Boling pt: Melting pt: Solubility: DMSO: 85 mg/mL (199.76 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Van Rompaey L, et al. J Immunol. 2013, 191(7), 3568-3577.

Note: This product is for in vitro research use only