


Description

Product Name	NVP 231
Brief Description	Inhibitors
Purification	98.00%
Target Name	ERK inhibitor
Calculated MW	431.55
Formulation	C25H25N3O2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:O=C(NC1=NC2=C(S1)C=C(NC(=O)C13CC4CC(CC(C4)C1)C3)C=C2)C1=CC=CC=C1Pathways:MAPKReceptor:CERKBoling

pt:Melting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Graf C, et al.Mol Pharmacol. 2008 Oct;74(4):925-32.

Note: This product is for in vitro research use only