

## Spatin-1

Catalog No: #S1937


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Spatin-1
Brief Description	Inhibitors
Purification	98.00%
Target Name	DUB inhibitor
Calculated MW	271.26
Formulation	C15H11F2N3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: FC1=CC=C(CNC2=NC=NC3=C2C=C(F)C=C3)C=C1  
Pathways: Ubiquitination Receptor: USP10; USP13  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: Faint Yellow Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Shao S, et al. Int J Oncol. 2014 May;44(5):1661-1668.

Note: This product is for in vitro research use only